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# Blood Pressure Monitor

*with Bluetooth and data management  
in a dedicated app*

## Operating Manual

BPM 1401  
adeVital  
Pressure

Thank you,

for deciding to purchase the **ADE Blood Pressure Monitor *adeVital BPM 1401 with App***. All products from the innovative *adeVital* program – body analysis scales, blood pressure monitor, and activity monitor – are equipped with a Bluetooth function for transmitting all data to the *adeVital-plus* app. The simple interface of the free downloadable app enables easy and clear monitoring of your progress and results. Health values are therefore available at any time and are always up-to-date.

## Table of Contents

<b>SAFETY AND USAGE INFORMATION</b> .....	<b>3</b>
<b>SAFETY INFORMATION</b> .....	<b>3</b>
USAGE INFORMATION .....	3
<b>FEATURES AND COMPONENTS</b> .....	<b>4</b>
FEATURES.....	4
THE LCD.....	5
COMPONENTS OF THE BLOOD PRESSURE MONITOR .....	6
<b>PACKAGING CONTENTS</b> .....	<b>7</b>
<b>START-UP</b> .....	<b>7</b>
INSERTING BATTERIES.....	7
EXTERNAL POWER SUPPLY.....	8
INSTALLING THE ADEVITAL PLUS APP AND CONNECTING.....	8
MONITORING FOR BLUETOOTH-ENABLED MONITORS .....	8
<b>PROGRAMMING AND CONNECTION/INPUT OF PERSONAL INFORMATION</b> .....	<b>9</b>
COUPLING OF BLOOD PRESSURE MONITOR AND ADEVITAL PLUS APP.....	9
<b>INDICATORS ON THE DISPLAY</b> .....	<b>9</b>
<b>COUPLING OF THE ADEVITAL PLUS APP WITH THE APPLE HEALTH APP</b> .....	<b>9</b>
<b>SETTINGS ON THE BLOOD PRESSURE MONITOR</b> .....	<b>10</b>
<b>STARTING MEASUREMENTS</b> .....	<b>11</b>
PUTTING ON THE CUFF .....	11
STARTING MEASUREMENTS .....	12
RETRIEVE STORED MEASUREMENT VALUES.....	12
MEASURING TIPS .....	13
<b>DATA TRANSMISSION</b> .....	<b>14</b>
<b>ERROR MESSAGES AND THEIR CORRECTION</b> .....	<b>14</b>
<b>ERROR MESSAGES</b> .....	<b>14</b>
DURING DATA TRANSMISSION.....	15
<b>MAINTENANCE AND CLEANING</b> .....	<b>15</b>
<b>YOUR BLOOD PRESSURE MONITOR AND ITS ENVIRONMENT</b> .....	<b>16</b>
<b>CE CONFORMITY DECLARATION</b> .....	<b>16</b>
<b>WARRANTY AND GUARANTEE</b> .....	<b>17</b>
<b>DISPOSAL</b> .....	<b>17</b>
<b>LEGAL DUTY TO INFORM ABOUT BATTERY DISPOSAL</b> .....	<b>17</b>
<b>DISPOSAL OF WASTE ELECTRONIC DEVICES</b> .....	<b>17</b>
<b>ABOUT BLOOD PRESSURE</b> .....	<b>17</b>
SYSTOLIC/DIASTOLIC BLOOD PRESSURE .....	17
STANDARD BLOOD PRESSURE CLASSIFICATION.....	18

IRREGULAR HEARTBEAT .....	18
<b>TECHNICAL DATA.....</b>	<b>19</b>
AVAILABLE ACCESSORIES .....	20
<b>EMC INFORMATION .....</b>	<b>20</b>

## SAFETY AND USAGE INFORMATION

### SAFETY INFORMATION

Warning signs and symbols are important for the correct and safe use of this product, and for protecting you and third parties from injury. The warning signs and symbols, which you will find on the labels and in this user manual, have the following meanings:

Symbol for: "Read the user manual!"

**Bluetooth**

Symbol for: "Bluetooth-enabled device"

Symbol for: "Type BF applied part"

**SN**

Symbol for: "Serial number"

Symbol for: "Direct current"

Symbol for: "Manufacturer"

For indoor use only

**F1**

T1A/250V  $\Phi$  3.6\*10CCC

Symbol for: "Protection Class II device"

Symbol for: "Complies with the requirements of European Directives"

Symbol for: "Environmental protection" – electrical products must not be disposed of with other household waste. Please recycle this device. Your local authorities or dealers will give you tips on recycling.

### USAGE INFORMATION

- The **ADE Blood Pressure Monitor adeVital BPM 1401 with App** enables blood pressure monitoring, pulse monitoring, and automatic storage of your measurement results. This blood pressure monitor has a lifetime of at least 3 years.
- The measurement results of the BPM1401 are equivalent to the results that a trained professional obtains with a cuff and stethoscope.
- This user manual contains important information relating to safety and maintenance, and it offers step-by-step instructions on the use of this product. Read this manual carefully prior to use.

## ATTENTION!

Please read this user manual carefully and attentively before using the device.

This device is designed for use exclusively by adults in the home environment.

This device is intended for the non-invasive measurement and monitoring of arterial blood pressure. It is not intended for use on extremities other than the upper arm, or for functions other than measuring blood pressure.

Do not confuse self-monitoring with self-diagnosis. With this device, you can monitor your blood pressure. Commence or terminate medical treatment only after consultation with your doctor. If you are taking medication, speak to your doctor about the appropriate time for measuring your blood pressure. Never change prescribed medications without the agreement of your doctor.

This device is not suitable for continuous blood pressure monitoring during medical emergencies or operations.

If the cuff pressure exceeds 40 kPa (300 mmHg), the cuff will be deflated automatically. Should the cuff fail to deflate when its pressure exceeds 40 kPa (300 mmHg), remove the cuff from your upper arm and press the appropriate button in order to stop the inflation.

Do not use the blood pressure monitor in conditions with strong electromagnetic fields.

The blood pressure monitor is not approved for categories AP or APG. It is not suitable for use in the presence of a flammable anaesthetic mixture with air (or oxygen, nitrous oxide).

Not suitable for children under 36 months of age, owing to the presence of small parts that may be swallowed.

Please use accessories and add-on parts that are recommended or authorised by the manufacturer, otherwise damage to the device may occur and the user may be put at risk.

Before using the blood pressure monitor, please ensure that the device functions safely and is in good working order.

## FEATURES AND COMPONENTS

### FEATURES

- User memory for 60 x 2 measurements on the device; unrestricted management of measurement results in the app
- Data transmission to the adeVital<sup>plus</sup> app via Bluetooth
- Precise measurement of blood pressure data (systolic and diastolic values) on the upper arm
- Measurement of pulse rate
- Indication of irregularities in the heartbeat (arrhythmia warning)
- Measurement during inflation
- Oscillometric measurement process, i.e. through the oscillation of the pulsing blood
- Classification of measured values according to pre-defined values compliant with the WHO\*
- Indication of date and time
- Automatic power-off
- Battery status indicator

THE LCD


SYMBOL	DESCRIPTION	EXPLANATION
	Data retrieval/Configuration button	Retrieve stored measurement values Set button for time setting
	User 1 button	Start measurement for User 1 and automatically transmit measurement results Programming button
	User 2 button	Start measurement for User 2 and automatically transmit measurement results Connection button for the adeVital <sup>plus</sup> app
<b>SYS</b>	Systolic blood pressure	Upper blood pressure value
<b>DIA</b>	Diastolic blood pressure	Lower blood pressure value
<b>Pul/min</b>	Pulse	Pulse beats/minute
	Pulse recognition	Displays your heartbeat, as long as it is recognised
	Data ready for transmission	Measurement data is stored in the device
	Data transmission	Data transmission in progress
	Average value	Average value of the last three measurements
	Bump memory	Bumps lead to imprecise measurements
	Batteries low	The batteries are low and have to be replaced
<b>mmHg</b> <b>kPa</b>	Unit	Unit of measurement of blood pressure
	Current time	Month: Day (Hour, Minute)
	Deflating	The cuff is being deflated
	Irregular heartbeat	Please read further information on page 16
	Bluetooth connection	Connects the BPM 1401 with the active app


## COMPONENTS OF THE BLOOD PRESSURE MONITOR

The following components are constituent parts of the device or are built into the device:


- ADE Blood Pressure Monitor adeVital BPM 1401, housing made of ABS plastic
- Soft, adjustable cuff: 22 – 42 cm
- Amplifier
- Air tube
- Pump
- Valve
- 128 mm x 50 mm blue LCD with white backlight
- Socket for connection to adaptor
- Battery compartment

## PACKAGING CONTENTS

– Blood pressure monitor  
BPM 1401

– 4 x AAA alkaline batteries 1.5 V

– Cuff  
(Type BF applied part)  
(AC2232-02)

– User manual

## START-UP

Take all of the components out of the packaging and remove the protective wrapping from around the components.


Keep the packaging for storage and/or possible return of the device.

This device is ***not*** intended for use by people (including children) with limited physical, sensory or mental capacities or lacking experience and knowledge, except when they are supervised by persons who are responsible for their safety or who are instructed in the use of the device. Children should not be allowed to play with the device.

## INSERTING BATTERIES

- Open the cover of the battery compartment on the rear side of the ADE Blood Pressure Monitor *adeVital* BPM1401
- Insert the batteries (4 x 1.5V AAA) into the battery compartment according to the polarity markings inside the compartment.
- Close the cover of the battery compartment.

### ATTENTION!

**REPLACE THE BATTERIES UNDER THE FOLLOWING CIRCUMSTANCES:**

appears on the LCD.

- The LCD becomes dimmer.
- When the blood pressure monitor is turned on, the LCD does not light up.
- Used batteries belong in hazardous waste. Do NOT dispose of these with other household waste. Please respect your local regulations and recycling advice when disposing of used batteries.

Battery capacity: 600 mAh. For three blood pressure measurements per day, 35 s per measurement, 20 s display of the measurement results, and 10 s data transmission. The current drawn during measurements totals 400 mA; during the display of measured data, 50mA; during data transmission, 50 mA (additional). The device, when switched off, requires 35  $\mu$ A.

Remove the batteries if the device is not going to be used for an extended period of time.


## EXTERNAL POWER SUPPLY

100 - 240 V~, 50 - 60 Hz, 400 mA (Only with mains adaptor UE08WCP-060100SPA! /Not Included!)

If you wish to stop mains operation, disconnect the cable of the mains adaptor from the blood pressure monitor.

Please note:

The socket for the mains adaptor is located on the right side of the device. So that you can easily remove the cable of the mains adaptor, do not place any objects on this side of the device.

## INSTALLING THE ADEVITAL PLUS APP AND CONNECTING

Using Bluetooth 4.0 technology, mobile devices with a Bluetooth function can receive your personal data. The adeVital<sup>plus</sup> app is available for free at the App Store and the Google Play Store. There, you can search for “adeVital plus” and install this app on your mobile device, or simply scan one of the QR codes appropriate to your device.

## MONITORING FOR BLUETOOTH-ENABLED MONITORS

- Display of daily data for all ascertained values
- Display of systolic and diastolic blood pressure values
- Pulse rate
- Generation of progress curves for various intervals
- Clear presentation and progress monitoring of all data and values as a diagram, graph, or table.
- Simple data transmission via integrated email function, Facebook, or Twitter.

Simply install the specially developed adeVital<sup>plus</sup> app and connect your scales to your mobile device. You can then make use of the comprehensive software in the adeVital<sup>plus</sup> app.

The adeVital<sup>plus</sup> app is available for iPod touch (5th Generation), iPhone 4S, iPhone 5, iPhone 5s, iPhone 6, iPhone 6 plus, iPad 3, iPad Air, iPad mini and later versions, as well as for Android devices from Version 3.0 up. Bluetooth 4.0 compatibility is required.

### ATTENTION!

BLUETOOTH MODULE NO.:	NRF8001
HF FREQUENCY RANGE:	2402 MHz to 2480 MHz
ACOUSTIC OUTPUT RANGE:	0 dBm
SUPPLY VOLTAGE:	3.3 V
TRANSMISSION RANGE:	10 m

## PROGRAMMING AND CONNECTION/INPUT OF PERSONAL INFORMATION

### COUPLING OF BLOOD PRESSURE MONITOR AND ADEVITAL PLUS APP

Open the installed adeVital<sup>plus</sup> app and record your personal information in all of the given fields. Under "Devices", select the large "+ADD DEVICE" symbol.

Do NOT switch the blood pressure monitor on! Press the "User 2" button for approx. five seconds until two alternately blinking zeroes appear on the display. The adeVital<sup>plus</sup> app will notify you as soon as your blood pressure monitor has been found. If no device has been found, press "update".


When the blood pressure monitor has been found, select the device in the adeVital<sup>plus</sup> app and couple it with the desired memory cell of the blood pressure monitor (1 or 2). Upon successful coupling, the [ ] symbol appears on the display of the blood pressure monitor.

Afterwards, confirm the input of your data.

You can set up the blood pressure monitor for 2 users, each with their own mobile device. For the second person who would like to use the blood pressure monitor together with the adeVital<sup>plus</sup> app, simply repeat the process described above and choose the free memory cell for this person.

### INDICATORS ON THE DISPLAY

– The 0 0 and 0 0 symbols appear alternately on the LCD of the blood pressure monitor and indicate that device coupling is in progress.

– Upon successful connection, the [ ] symbol appears on the LCD.


– In the event of a failed connection, "E1" appears on the LCD.

### COUPLING OF THE ADEVITAL PLUS APP WITH THE APPLE HEALTH APP

If you use the Apple Health app, you can connect the app to the adeVital<sup>plus</sup> app and then analyse your values in the Apple Health app, too.

To couple both apps, open the "Management" menu in the adeVital<sup>plus</sup> app and go to "Settings" => "Health-Kit". In the "Vital Signs" category, activate the data (blood pressure, heart rate) that you would also like to analyse in the Apple Health app. Then open the Apple Health app and, in the "Source" menu, select the "adeVital plus" app. Here, in the "Vital Signs" category, activate the data according to your selection in the adeVital<sup>plus</sup> app.

After coupling the adeVital<sup>plus</sup> app and the Apple Health app, you can analyse all subsequently measured values in the adeVital<sup>plus</sup> app as well as in the Apple Health app.

## SETTINGS ON THE BLOOD PRESSURE MONITOR

Please set the time before using the blood pressure monitor for the first time, so that all recorded measurement data is furnished with the correct timestamp. (Default value – Date: 2012-01-01; Time: 00:00; Unit of measurement: mm Hg)

In settings mode, the blood pressure monitor automatically deactivates 60 seconds after the last input.

1. When the blood pressure monitor is switched off, press and hold the "User 1" (SET) button in order to set the time and date.


2. The hours and minutes indicator will appear on the display with a blinking hour value.  
Set the hour value by pressing the button; this will change from **AM** to **PM** after 12 hours have elapsed.

3. By pressing the **SET** button, change to minutes and set these using the button. By holding down the button, you can cycle more quickly through the numbers.

4. Repeat steps 2 and 3 in order to confirm the month (M) and day (D) settings.

5. After pressing the **SET** button (to confirm the day), the year value will blink.

6. Hold down the button until the correct year value appears. The values appear on a loop; after [2052] the numerical sequence begins again with [2012]. The confirm the [YEAR], press the **SET** button and switch to [Unit of measurement].

7. Repeat steps 2 and 3 in order to confirm the [UNIT OF MEASUREMENT] setting.

[mmHg]

[kPa]


8. After completing and saving the inputs, the word **[DONE]** appears on the display and the device deactivates.


## STARTING MEASUREMENTS

### PUTTING ON THE CUFF

1. Stick the plug of the air tube on the cuff into the tube connection socket on the left side of the blood pressure monitor.


2. Put the cuff on your upper arm (see image on right).

3. Close the cuff. Ensure that the cuff is located 2 to 3 cm above your elbow.

When the cuff has been closed around your upper arm, you should only be able to slide one finger between the cuff and your arm.

4. Correct posture:

Only wear close-fitting clothing or expose your arm while measuring.

Sit in a comfortable and relaxed position on an appropriate chair.

The middle of the cuff should sit level with your heart.

Legs relaxed, feet pointing outwards.

The palm of your hand should point upwards.

## STARTING MEASUREMENTS

After correctly position the cuff, press the "User 1" ("User 2") button in order to switch on the blood pressure monitor. The measurement process will then be carried out automatically.

1. LCD

2. Automatic zeroing

3. Automatic inflation and measurement.

4. Display and saving of results, which are automatically transmitted to the adeVital<sup>plus</sup> app. The «» symbol blinks.

5. Press "User 1" ("User 2") in order to switch off the blood pressure monitor. Alternatively, the blood pressure monitor will switch itself off within 1 minute of the last actuation process.

A. Once the entire measurement has been completed, press a different button (button for User 1 or 2); the blood pressure monitor will then start to take a new measurement.

B. A maximum of 60 measurement results can be recorded for User 1 and User 2 on the device, and unlimited measurement results in the adeVital<sup>plus</sup> app.

## RETRIEVE STORED MEASUREMENT VALUES

1. When the blood pressure monitor is switched off, press the "Search" button in order to retrieve the saved measurement values. The blood pressure monitor displays the average value of the last three measurements.

2. Press the "Search" button again in order to scroll through the records. Up to 60 records are saved for each user ID.

The date and time of the measurements will be alternately displayed.

3. Press the "User 1" and "User 2" buttons in order to switch between users.

NO RECORD FOR USER 2

4. Press "User 1" ("User 2") in order to switch off the blood pressure monitor. Alternatively, the blood pressure monitor will switch itself off within 1 minute of the last actuation process.

### MEASURING TIPS

UNDER THE FOLLOWING CIRCUMSTANCES, IMPRECISE MEASUREMENTS MAY ARISE AND SHOULD BE AVOIDED:

Within 1 hour after eating or drinking

Measuring immediately after consuming tea or coffee, or after smoking.

Within 20 minutes after a bath.

When talking or when moving your fingers

In a very cold environment

When your bladder is overfull

## DATA TRANSMISSION

Upon successful coupling of the blood pressure monitor *adeVital* BPM1401 and your mobile device, the measurement data will be automatically transmitted via Bluetooth to your mobile device, when the app is open.

The  symbol disappears following a successful data transmission. You can then check and analyse your personal health data on your mobile device.

In the event that the data transmission fails, the  symbol will be displayed again. The pending measurement data will be stored temporarily in the blood pressure monitor and transmitted to your mobile device upon completion of the next measurement.

### ATTENTION!

- In the vicinity of devices that bear the  symbol, interference is possible. The blood pressure monitor can interfere with electrical devices in the vicinity.
- For data transmission, this product must be paired with a Bluetooth device at 2.4 GHz.

#### How can potential interference be reduced?

1. The distance between the blood pressure monitor and the Bluetooth device should be as small as possible, between 1 m and 10 m. No obstacles should impair the Bluetooth connection between the blood pressure monitor and the Bluetooth device.
2. To prevent interference, other electronic devices (in particular, those with Bluetooth-transmission/transmitters) should be kept at a distance of at least 1 m from the blood pressure monitor.

## ERROR MESSAGES AND THEIR CORRECTION

This section contains a list with error messages and frequently asked questions relating to problems that can arise with your blood pressure monitor. If the device does not work as expected, please check this list before you send the device in for repair.

### ERROR MESSAGES

Problem	Symptom	Please check	Remedy
No current	The display is weak or does not light up.	Batteries have been depleted.	Insert new batteries.
		Batteries have not been correctly inserted.	Insert batteries correctly.
Batteries low	 on the display	Batteries are low.	Insert new batteries.

Error message	Please check	Remedy
E 1 appears	The cuff is loose.	Tighten the cuff and then take a new measurement.
E 2 appears	The cuff is too tight.	Tighten the cuff and then take a new measurement.
E 3 appears		Relax for a moment and then take a new measurement.
E 10 or E 11 appears		Movements can affect the measurement. Relax for a moment and then take a new measurement.
E 20 or E 21 appears		Relax for a moment and then take a new measurement.
Eexx appears on the display	A calibration error has occurred.	Repeat the measurement. If the problem persists, please contact the dealer or our customer service division. You will find contact information and return shipping instructions in the warranty section.

#### DURING DATA TRANSMISSION...

Error	Description	Solution
Data transmission failed.	Bluetooth is OFF (deactivated).	Switch Bluetooth on via "Setting >> General >> Bluetooth".
	The app is OFF (deactivated).	Tap the symbol of your app in order to start it.
	Outside the maximum Bluetooth range.	Place your iPhone closer to the blood pressure monitor.

#### MAINTENANCE AND CLEANING


- Prior to use
- Take the entire device out of the storage bag. Use a soft cloth to remove dirt from the blood pressure monitor. Clean the cuff with alcohol before attaching it.
- After use
- Wipe down the device using a soft cloth. Disinfect the cuff with some alcohol. Then place the entire device into the storage bag.


Please note that the right to exchange installed or attached parts is reserved **ONLY FOR SERVICE STAFF OF THE MANUFACTURER.**

## YOUR BLOOD PRESSURE MONITOR AND ITS ENVIRONMENT


Please bear in mind the following advice for the use of the blood pressure monitor:


Keep the device in a dry location and avoid direct sunlight


Do not immerse the device in water

Avoid shocks and impacts

Avoid dusty environments and temperature fluctuations

Remove dirt using a lightly moistened cloth

Do not wash the cuff

- Proper care significantly improves the service life of your blood pressure monitor.
- Do not use any aggressive cleaning agents, or abrasive or other chemicals to remove dirt from the device. This can lead to colour changes, malfunctions, or damage.
- The device must not be disassembled or modified. If you have problems, please contact your dealer or our customer service division.

## CE CONFORMITY DECLARATION

- This product has been manufactured in accordance with the harmonised European standards. They conform to the regulations of the EC-directives listed below:
- EMC Directive 2004/108/EC
- R&TTE Directive 99/5/EC
- MDD Directive 93/42/EWG

## WARRANTY AND GUARANTEE

Waagen-Schmitt guarantees the cost-free correction of defects due to material or manufacturing faults through repair or exchange (excluding batteries and parts subject to wear) for 3 years from the date of purchase. Please get your dealer to fill out and stamp the warranty section at the time of purchase. In the case of a warranty claim, please return the blood pressure monitor with the warranty section specifying the reason for complaint to your dealer.

**Waagen-Schmitt GmbH**

Hammer Steindamm 27–29, 22089 Hamburg

## DISPOSAL

### LEGAL DUTY TO INFORM ABOUT BATTERY DISPOSAL

Batteries do not belong in household waste. As a consumer you are legally obligated to return batteries. You can hand in your old batteries at the public collection point in your community or anywhere where batteries of the same kind are sold.

#### Please note:

These symbols can be found on batteries:

Li = Battery contains lithium

Al = Battery contains alkali

Mn = Battery contains manganese

CR (Li); A (Al, Mn); AA (Al, Mn); AAA (Al, Mn)

### DISPOSAL OF WASTE ELECTRONIC DEVICES


This symbol on the product or its packaging indicates that the product is not to be treated like normal household waste, but must be handed over at a collection point for the recycling of electrical and electronic devices.

Further information can be obtained in your community, at communal disposal services or the business where you bought the product.

## ABOUT BLOOD PRESSURE

### SYSTOLIC/DIASTOLIC BLOOD PRESSURE

When the chambers of the heart contract and blood is pumped out of the heart, the blood pressure attains its maximum value in the cycle. This value is called the systolic pressure. When the chambers of the heart relax, the blood pressure attains its minimum value in the cycle. This value is called the diastolic pressure.


## STANDARD BLOOD PRESSURE CLASSIFICATION

The blood pressure classification published in 1999 by the World Health Organisation (WHO) and the International Society of Hypertension (ISH) appears as follows:

### ATTENTION!

Only a doctor can tell you your normal blood pressure range. Please visit a doctor if your blood pressure lies outside the range. Please keep in mind that only a doctor can tell you whether your blood pressure has reached a dangerous level.

Value Blood pressure (mm Hg)	Low	Normal	High-normal	Mild	Moderate	Severe
SYS	<100	101-129	130-139	140-159	160-179	≥180
DIA	<60	61-84	85-89	90-99	100-109	≥110

## IRREGULAR HEARTBEAT

This blood pressure monitor is equipped with an intelligent IHB (*Irregular Heartbeat Detector*) for recognising irregular heartbeats. With every measurement, the IHB charts your heartbeat intervals and ascertains the standard deviation. If the calculated value is greater than or equal to 15, the IHB symbol appears on the display together with the measurement result.

### ATTENTION!

The IHB symbol indicates that during the measurement a pulse irregularity was detected in conjunction with an irregular heartbeat. This is normally NOT a reason for concern. However, if the symbol is frequently displayed, we recommend that you seek medical advice. Please bear in mind that this device is not a substitute for any investigation of the circulatory system, but can detect pulse irregularities at an early stage.

## DAILY VARIATIONS IN BLOOD PRESSURE

- Individual blood pressure fluctuates during the day. It also depends on how tightly you put the cuff on and in what position you take the measurements. You should therefore ensure that you take measurements under similar conditions.
- The variations in blood pressure are larger when a person is on medication.
- Wait at least 10 minutes before the next measurement.

## DIFFERENT BLOOD PRESSURE VALUES IN HOSPITAL OR AT THE DOCTOR AND AT HOME.

Blood pressure varies over the course of the day. Reasons for this are weather, emotions, movement, etc., and especially the “white coats” in hospital or at the doctor, which lead to higher blood pressure values than at home.

## WHEN TAKING YOUR BLOOD PRESSURE AT HOME, YOU MUST ENSURE THE FOLLOWING:

- The cuff must be correctly put on.
- The cuff must not be too tight or too loose.
- The cuff must not be put on the wrist.
- You must not be anxious when the cuff is being inflated.
- You should take 2 to 3 deep breaths before the measurement begins.

Tip: Take a few minutes to prepare for the measurement, until you are calm.

## DIFFERENT BLOOD PRESSURE VALUES ON BOTH OF MY ARMS

Blood pressure can be measured on both arms, but the measurements results will be different on each arm. You should therefore measure your blood pressure on the same arm each time.

## TECHNICAL DATA

### Power supply:

Battery operation:	6 V (4 x AAA alkaline batteries)
Mains operation:	100 - 240 V~, 50 - 60 Hz, 400 mA (only with mains adaptor UE08WCP-060100SPA!)

### Display mode:

Blue LCD with white backlight  
Display surface: 128 mm (L) x 50 mm (W)

### Measurement mode:

Oscillographic test method

### Measuring range:

Blood pressure: 0 - 40 kPa (0 ~ 300 mmHg)  
Pulse value: (40 - 199) pulse beats/minute

### Accuracy of blood pressure:

At 5 °C – 40 °C:	within ± 0.4 kPa (3 mm Hg)
At 0 °C – 45 °C (but excluding 5 °C – 40 °C):	within ± 0.7 kPa (5 mm Hg); Pulse value: ±5%

### Operating conditions:

Temperature:	5 °C – 40 °C
Relative humidity:	≤85%
Air pressure:	86 – 106 kPa

### Storage and Transport Conditions:

Temperature:	-20 °C – +60 °C
Relative humidity	10% - 93%
Air pressure:	50 – 106 kPa

### Measurable circumference of upper arm:

approx. 22 – 42 cm

### Weight:

approx. 385 g (without batteries)

### External dimensions:

approx. 120 x 160 x 69 mm

### Accessories:

4 x AAA alkaline batteries, user manual

### Mode of operation:

Continuous operation

### Protection class:

Type BF applied part

### Device classification:

Battery operation: ME device with internal power supply via  
mains adaptor: Class II ME device

### IP classification:

IP22

### Software Version:

V01

## AVAILABLE ACCESSORIES

<b>Mains adaptor</b> <b>Adaptor type:</b> UE08WCP-060100SPA <b>Input:</b> 100 - 240 V, 50 - 60 Hz, 400 mA <b>Output:</b> 6 V $\overline{\text{---}}$ 1 A	
---	---

Please visit our shop: [www.waagen-schmitt-shop.de/adevital/](http://www.waagen-schmitt-shop.de/adevital/)

## EMC INFORMATION

MANUFACTURER DECLARATION – ELECTROMAGNETIC EMISSIONS (IEC 60601-1-2)		
The Blood Pressure Monitor BPM1401 is intended for operation in an electromagnetic environment, as specified below. The user of the Blood Pressure Monitor BPM1401 should ensure that it is operated in an environment of this kind.		
Emitted interference measurements	Compliance	Electromagnetic environment – guidelines
HF emissions according to CISPR 11	Group 1	The Blood Pressure Monitor BPM 1401 uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and it is unlikely that adjacent electronic devices will be subject to interference.
HF emissions according to CISPR 11	Class B	The Blood Pressure Monitor BPM1401 is suitable for use in all establishments, including the home as well as establishments that are directly connected to the public low-voltage mains network which supply private households.
Emission of harmonics according to IEC 61000-3-2	Not applicable	
Emission of voltage fluctuations/flicker according to IEC 61000-3-3	Not applicable	

MANUFACTURER DECLARATION – ELECTROMAGNETIC IMMUNITY (IEC 60601-1-2)			
The Blood Pressure Monitor BPM1401 is intended for operation in an electromagnetic environment, as specified below. The customer or user of the Blood Pressure Monitor BPM1401 should ensure that it is operated in an environment of this kind.			
Immunity Tests	IEC 60601 Test Level	Compliance Level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) according to IEC 61000-4-2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be made of wood or concrete, or covered with ceramic tiles. If the floor is covered with synthetic material, the relative air humidity must be at least 30%.
Electrical fast transient / burst test according to IEC 61000-4-4	±2 kV for power cable	±2 kV for power cable	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Surges according to IEC 61000-4-5	±1 kV cable-cable	±1 kV normal-mode voltage	The quality of the supply voltage should correspond to that of a typical business or hospital environment.

Voltage drops, short-term interruptions, and fluctuations of the supply voltage according to IEC 61000-4-11	<5% $U_T$ (>95% drop in $U_T$ ) for 1/2 period	<5% $U_T$ (>95% drop in $U_T$ ) for 1/2 period	The quality of the supply voltage should correspond to that of a typical business or hospital environment. If the user of the Blood Pressure Monitor BPM1401 requires continued functioning even when interruptions of the energy supply occur, it is recommended that the Blood Pressure Monitor BPM1401 be supplied from an interruption-free power supply or a battery.
	40% $U_T$ (60% drop in $U_T$ ) for 5 cycles	40% $U_T$ (60% drop in $U_T$ ) for 5 cycles	
	70% $U_T$ (30% drop in $U_T$ ) for 25 cycles	70% $U_T$ (30% drop in $U_T$ ) for 25 cycles	
	<5% $U_T$ (>95% drop in $U_T$ ) for 5 s	<5% $U_T$ (>95% drop in $U_T$ ) for 5 s	
Magnetic field at the supply frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	The line-frequency magnetic fields should correspond to the characteristics of a typical installation site in a commercial or clinical environment.
<b>NOTE:</b> $U_T$ is the alternating mains voltage prior to the application of the testing level.			

<b>MANUFACTURER DECLARATION – ELECTROMAGNETIC IMMUNITY (IEC 60601-1-2)</b>			
The Blood Pressure Monitor BPM1401 is intended for operation in an electromagnetic environment, as specified below. The user of the Blood Pressure Monitor BPM1401 should ensure that it is operated in an environment of this kind.			
Immunity Tests	IEC 60601 Test Level	Compliance Level	Electromagnetic environment – guidelines
Grid-bound HF interference immunity according to IEC 61000-4-6  Radiated HF interference immunity according to IEC 61000-4-3	3 Vrms 150 kHz up to 80 MHz  3 V/m 80 MHz up to 2.5 GHz	3 Vrms  3 V/m	<p>Portable and mobile RF devices should not be used in closer proximity to the Blood Pressure Monitor BPM1401 (including cables) than the recommended protective distance, which is calculated according to the appropriate equation for the transmission frequency.</p> <p><b>Recommended protective distance:</b></p> $d = 1.167 \sqrt{P}$ $d = 1.167 \sqrt{P} \text{ for 80 MHz up to 800 MHz}$ $d = 2.333 \sqrt{P} \text{ for 800 MHz up to 2.5 GHz}$ <p>where P is the nominal power of the transmitter in watts (W) according to the manufacturer of the transmitter, and d is the recommended protective distance in metres (m).</p> <p>The field strength of stationary RF transmitters should, according to an examination on site<sup>b</sup>, be lower than the compliance level.<sup>b</sup></p> <p>In the vicinity of devices that bear the following pictorial symbol, interference is possible:</p> 

COMMENT 1	At 80 MHz and 800 MHz, the higher frequency range applies.
COMMENT 2	These guidelines may not be applicable in all cases. The proliferation of electromagnetic fields is affected by the absorptions and reflections of buildings, objects, and people.
<sup>a</sup>	<b>The field strengths of stationary transmitters, such as base stations of cordless telephones and mobile radios, amateur radio stations, AM and FM radio and television transmitters, cannot be theoretically predetermined with precision. In order to determine the electromagnetic environment in regard to stationary transmitters, a study of the electromagnetic phenomena of the site must be considered. If the measured field strengths on the site where the Blood Pressure Monitor BPM1401 is to be used exceed the compliance levels specified above, then the Blood Pressure Monitor BPM1401 should be studied in order to verify that it is functioning as intended. If unusual performance indicators are observed, additional measures may need to be taken, such as a modified orientation of or another location for the Blood Pressure Monitor BPM1401.</b>
<sup>b</sup>	<b>Over the frequency range from 150 kHz to 80 MHz, the field strength should be lower than 3 V/m.</b>

Recommended protective distances between portable and mobile HF telecommunication devices and the ME device or ME system – for ME devices or ME systems that are not life-sustaining.

<b>RECOMMENDED PROTECTIVE DISTANCES BETWEEN PORTABLE AND MOBILE HF TELECOMMUNICATION DEVICES AND THE BLOOD PRESSURE MONITOR BPM 1401.</b>			
The Blood Pressure Monitor BPM1401 is intended for operation in an electromagnetic environment in which the HF disturbances are controlled. The customer or user of the Blood Pressure Monitor BPM1401 can help to prevent electromagnetic interference by complying with the minimum distance between portable and mobile HF telecommunication devices (transmitters) and the Blood Pressure Monitor BPM1401 – depending on the output performance of the communication device, as specified below.			
Nominal power of the transmitter (W)	Protective distance, dependent on transmission frequency (m)		
	150 kHz to 80 MHz $d = 1.167 \sqrt{P}$	80 MHz to 800 MHz $d = 1.167 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.333 \sqrt{P}$
0.01	0.167	0.167	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.690	3.690	7.338
100	11.67	11.67	23.33
For transmitters whose maximum nominal power is not specified in the table above, the recommended protective distance $d$ in metres (m) can be determined by using the equation that pertains to a given column, where $P$ is the maximum nominal power of the transmitter in watts (W) according to the manufacturer of the transmitter.			
COMMENT 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
COMMENT 2	These guidelines may not be applicable in all cases. The proliferation of electromagnetic fields is affected by the absorptions and reflections of buildings, objects, and people.		

**Warranty section – 3 years on the blood pressure monitor** (excluding batteries and parts subject to wear)

• Sender


• Reason for complaint

<hr/>	<hr/>
<hr/>	<hr/>
	<hr/>
	<hr/>

• Date of Purchase \_\_\_\_\_

Please send this section incl. proof of purchase back to your specialist retailer.

Operating Manual\_BPM 1401\_EN\_130527\_REV003

*Collect, manage, and analyse health data: It's simple and clear using the new adeVital product series with Bluetooth functionality and the intelligent adeVital<sup>plus</sup> app! That's the nature of modern health!*

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