Bedienungsanleitung
Operating manual
Manuel d'utilisation
Manual de instrucciones

ADE

Klassische Säulenwaage mit Längenmessung
Classic beam scale with height rod
Balance classique à fléau avec toise
Báscula mecánica clásica con tallímetro

M318800

M110800-210630-Rev004-UM-de-e-es

Inhalt – Content – Sommaire – Contenido

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www.ade-germany.de

General

This symbol means: "Consult operating instructions". To achieve accurate results, read the instructions for use carefully before use and follow all instructions provided therein. Keep the instructions for use for future reference.

Your ADE Classic Physician Beam Scale with height rod is a precision instrument and must be treated accordingly. Therefore, kindly read and follow the operating manual in all respects.

Check carton contents:

- Column with head part
- Base/platform
- Height rod
- Pair of wheels (incl. 2 cross-slotted crews)
- Levelling feet 4 pieces
- Mounting material: (6 x cross-slotted screw, 2 x hexagonal screw, 1 wrench)
- Operating manual

Safety

- Handle your scale with care - it is a precision instrument.
- Use the scale on a firm and level surface
- Operate the unit only within the permissible ambient conditions.
- Do not operate the unit in hazardous or unstable environments.
- Clean with a dry cloth, do not use harsh cleaning agents.
- Do not immerse the appliance in water or other liquids.

Assembling


Fig. 1 - ADE M318800

The assembling of the scale is easily done as all parts and components are pre-assembled during manufacturing:

Fixing of levelling feet


Fig. 2 Levelling Feet

The levelling feet are to be screwed into the holes underneath the weighing platform.

Assembling of the pillar


Fig. 3 - Pillar assembling

1. Insert the pillar upright into the bracket
2. Use the screws and washers to fix the pillar in a vertical position

Connection of the steelyard rod to the long lever


Fig. 4 - Connection of the steelyard


Fig.4-1

Fig. 4-2


Fig. 4-3

1. Place the pillar on a table (Fig. 4).
2. Disassemble the transport wires of the levers.
3. Insert the wrench to the small hole in front of the steelyard rod and pull the steelyard rod hook out.
4. Push the long lever up and hook the bearing of the steelyard rod over the pivot of the long lever (Fig. 4-2).
5. Take the wrench off (Fig. 4-3).
6. Place the scale gently on the ground.


The connecting rod must move freely.

Height rod Installation

Height rod installation

Fig 5

Back side hole of the height rod

Fig 5-1


Fig 5-2


Fig 5-3

1. Place the height measure with the holes at the back side onto the bolts at the front side of the pillar, adjust the height rod downwards to the right position (Fig. 5-1).
2. Use the wrench to tighten the two hex-head screws. Do not over tighten the screws (Fig. 5-2).
3. Take the wrench away after tightening of the screws (Fig. 5-3).

Installation of castors – moving wheels



Fig 6-1

Moving wheel installation

Fig 6

1. Place the wheel-set with the two holes towards the frame (platform) of the scale. The holes must be in an upside position.
2. Use the screws and washers to fix the wheel-set to the frame (platform) and adjust the wheel-set at a level position, then tighten the screws.
3. When moving the scale first tilt the scale on the wheels and make sure, that the platform is well off the ground.

Correct weighing and measuring

Zero adjustment:

Before the first use of the scale, please adjust the zero point carefully to insure correct weighing results.

- First step on and off the scale for several times to make sure that all parts come back into the right position after transport.
- For zero setting keep the platform empty.
- Set both sliding weights on the beam (head part) into zero position.
- The beam has to move freely equilibrating very close to level.
- Turn the zero adjustment screw (Fig. 3) slowly to the right respectively to the left to bring the beam into the exact level position.
- The scale is now ready to be used.

A periodically zero setting according to the previous instructions should be done. This applies especially if the scale is used in a different place and on a different surface.

Weighing

- Step on the scale.
- Slide the larger counterweight to the estimated weight.
- Now use the small sliding weight to reach the correct result.
- Wait until the beam comes close to a standstill in even position.
- Even if the head is still moving you can read the correct result +/- 100 g.
- Total weight is the addition of both results (smaller and larger counterweight).

Length Measuring

- Place measuring tongue into a horizontal position and pull the height rod up.
- Place the person on the scale and place the measuring tongue on top of the head.
- For the person shorter than 114,7 cm push the red button on top of the measure and lower the measuring tongue until it reaches the person's head.
- The reading is done either at the red line or between the two white triangles.
- If the measure is not in use, place the measuring tongue in a vertical position and push the height measure down.

Cleaning and maintenance

Cleaning

For cleaning just use a damp cloth or common disinfection fluids.
Do not use any aggressive cleaning fluids.

Warranty – Liability

ADE will repair or replace this product within a period of 2 year from date of delivery applying to defects occurred due to poor material or workmanship (presentation of invoice is necessary). All moveable parts – batteries, cables, mains units, rechargeable batteries etc. – are excluded. Warranty does not include natural wear, damages caused by incorrect or negligent handling, overuse or unsuitable or improper storage. The cost of transport will be borne by the customer should the equipment be located anywhere other than the customer's premises. In the event of transport damage, claims under warranty can be accepted subject to using the complete original packaging and the equipment secured the way it was packed and secured originally. All the packaging should therefore be retained.

Any warranty claim will be rejected if the equipment is opened by persons not authorized expressly by ADE to do so.

Customers abroad are kindly requested to contact their local dealer for any warranty claim, customer service and spare parts.

Technical Data

	M318800	M318800-02
• Max. Capacity:	200 kg	200 kg / 440 lb
• Graduation:	100 g	100 g / 4 oz
• Platform size:	375 x 275 mm	
• Dimensions:	272 x 1467 x 530 mm	

Height rod:

- Measuring range: 600 – 2100 mm
- Graduation: 1 mm
- Total Weight: approx. 18 kg

CE-conformity

Declaration of conformity by the manufacturer

These scales have been manufactured in accordance with the harmonised European standards. They conform to the regulations of the EC-directives listed below:

- Machinery Directive 2006/42/EC
in the valid version.

This declaration loses its validity if modifications are made to the scale without our approval.

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ADE Germany GmbH