

# Electronic wheel chair and chair scales

## M500020 Electronic wheel chair scale

**ADE**



Low level platform scale with one ramp for convenient access with wheel chairs. A cable remote display shows the weight clearly and can be mounted to the wall. Two castors and one handle permit a quick transport. The pre tare function will allow you to deduct a memorized weight just by pressing one button. Such the net weight of a patient can be determined. A hold and tare function will help to get a quick and precise reading of the patient's weight. Powered by mains and rechargeable batteries.

### Technical specification:

Capacity:	300 kg
Graduation:	100 g
Dimensions:	915 x 800 x 44 mm
Weight:	22 kg
Platform size:	800 x 800 mm
Functions:	Automatic zero setting, tare, pre-tare, low battery indication, auto witch off, hold, BMI calculation
Option:	MZ30010 Second ramp