



This extra-large platform scale will take generally used stretchers of various sizes and allows to weight patients lying on a stretcher during transport. The use as wheel chair scale is a further option. The pre tare function will allow you to deduct a memorized weight just by pressing one button. Such the net weight of a patient can be determined.

Technical specification:

Capacity:	300 kg
Graduation:	100 g
Dimensions:	1515 x 800 x 65 mm
Weight:	28 kg
Platform size:	1400 x 800 mm
Functions:	Automatic zero setting, tare, pre-tare, low battery indication, auto witch off, hold, BMI, RS232-Interface