

Electronic wheel chair and chair scales

ADE

M500020-02 Electronic wheel chair scale with foldable column



Low level platform scale with one ramp for convenient access with wheel chairs and foldable display holder column. Two castors and one handle permit a quick transport. The pre tare function will allow you to deduct a memorized weight just by pressing one button. Such the net weight of a patient can be determined. A hold and tare function will help to get a quick and precise reading of the patient's weight. Powered by mains and rechargeable batteries.

Technical specification:

Capacity:	300 kg
Graduation:	100 g
Dimensions:	915 x 800 x 44 mm
Weight:	22 kg
Platform size:	800 x 800 mm
Functions:	Automatic zero setting, tare, pre-tare, low battery indication, auto witch off, hold, BMI calculation
Option:	MZ30010 Second ramp