

Bariatric, bed and lifter scales

0.0 kg



For the use of bedridden patients in intensive care units and an essential must in renal kidney units during dialysis treatment. The bed can be easily pushed on the 4 weighing shoes with the help of the integrated ramps. The inbuilt rechargeable battery will store the weight even if an unwanted power failure occurs. The instrument trolley will be used for storing the cells and for transport of the complete weighing scale.

Technical specification:

Capacity:	500kg total capacity patient weight 300kg
Graduation:	100 g
Dimensions:	600 x 550 x 920 mm
Weight:	33 kg
Functions:	Automatic zero setting, tare, pre-tare and memory, low battery indication, auto switch off, hold, BMI, RS232-Interface, including mobile stand

