## Electronic baby scale M101000-01 CE class III approved







Electronic baby and toddler weighing scale for clinics and hospital use. The removable soft curved baby tray gives comfortable space to babies and when removed allows toddlers to be weighed in a standing up position. The weight of infants is measured precisely in 5 g steps. Operated by standard 6 x AA batteries and mains. Optionally is a carrying bag MZ10060 available.

Technical specification:	
Capacity:	20 kg
Graduation:	5g<10kg>10g
Functions:	Tare, hold, automatic zero setting, low
	battery indication, automatic switch off
Dimensions:	560 x 350 x 130 mm
Dimensions t	ray: 560 x 280 x 95 mm
Weight appro	ox.: 3,4 kg

