

Electronic personal scale  
M320000-02 CE class III approved

**ADE**



Practical floor scale with a slightly upraised column for heavy duty professional use. Clear large LCD display with easy tap-on-start and auto-hold function for quick and precise weighing. The low level slip resistant platform offers a wide and comfortable space for the patients. Straight forward BMI, Auto-Hold and mother-child operation gives fast results of the patients' nutritional standard. The turning display allows easy reading from patient and user side. The scale operates by mains and batteries.

**Technical specification:**

Functions: Tare, hold, automatic zero setting, turning display, low battery indication, auto switch off, BMI calculation, kg, mother & child function

Capacity & graduation: 250kg by 50g / 100g

Dimensions & weight: 450x350 x920mm 7kg

