
Electronic column weighing scale M320000-01 CE class III approved







The perfect choice for use in all medical fields. Heavy duty electronic column scale with an extra wide, flat and slip resistant platform. The transport castors ensure easy movability. The BMI, hold and mother-child functions provide fast and clear results of the patients' nutritional condition. The turnable display allows easy reading from patient and user side. The scale operates by mains and batteries, included as standard. The scale setup is fast and easy: supplied in just two parts and assembled by one set of special screws.

Technical specification:

Functions:	Tare, hold, automatic zero setting, turning display, low battery indication, auto switch off,
	BMI calculation, kg, mother & child function
Capacity & graduation:	50g <50kg> 100g
Dimensions & weight:	350 x 450 x 920 mm; 7,7 kg
Optionally	MZ10023-1 telescopic length measure

