

Electronic wheel chair scales

M500020-03 CE class III and MDD approved

ADE



Low level platform scale with one ramp for convenient access and a foldable heavy-duty handrail. Two castors and one handle permit a quick transport. Powered by mains and rechargeable batteries. Allowing the use of the scale at any place. A second ramp is available as option.

Technical specification:

Capacity:	300kg total capacity patient
Graduation:	100 g
Dimensions:	915 x 800 x 900mm
Weighing platform:	800 x 800 mm
Weight;	34kg
Power supply:	By mains and rechargeable batteries
Functions:	Automatic zero setting, tare, pre-tare, low battery indication, auto witch off, hold, BMI, RS232-Interface
Option:	MZ30010 second ramp