

Electronic chair weighing scale  
M403020 CE class III approved

**ADE**

## Wheelchair and chair scales

0.0 kg



Due to its 4 steering wheels of which are 2 fitted with brakes, this chair scale can be moved into the smallest corner without any effort. The foldable footrests allow a trouble free access for the patient. The padded arm rests and the ergonomically seating position make the scale comfortable to use.

### Technical specification:

Capacity:	250 kg
Graduation:	100 g
Dimensions:	600 x 850 x 960 mm
Weight;	28 kg
Seat size:	380 x 450 mm / seat height 490 mm
Functions:	Automatic zero setting, tare, pre-tare, low battery indication, auto witch off, hold, BMI, RS232-Interface