

Electronic chair weighing scale
M400020 CE class III approved

ADE



Comfortable and practical to use the chair scale. Swivelling foot and arm rests allow a comfortable access in and out of the chair. Four wheels make the scale easily manoeuvrable and safe due to two brakes at the back wheels. Tare and the hold function help to handle all patients. Powered by rechargeable batteries or mains for the use at any place.

Technical specification:

| | |
|-------------|--|
| Capacity: | 250 kg |
| Graduation: | 100 g |
| Seat size: | 420 x 480 mm / seat height 530 mm |
| Functions: | Automatic zero setting, tare, pre-tare and memory, low battery indication, auto switch off, hold, BMI, RS232-Interface |