



Electronic baby hanging scale - small and practical device made of coated aluminum, with integrated length measuring tape and a fixed hook made of galvanized metal for weighing slings. The large folding handle made of blasted stainless steel provides an extra secure hold and allows newborns to be weighed quickly and safely. The scale weighs babies up to 20 kg with an accuracy of 20 g. The LCD with blue backlight enables clear reading even in low light conditions.

### Technical specification:

|              |  |
|--------------|--|
| Capacity:    | 20 kg  |
| Graduation:  | 20 g   |
| Functions:   | Tara, unit switch, automatic zero setting, low battery indication, auto switch off, auto-hold function |
| Dimensions.  | 55 x 22 x 100 mm   |
| accessories. | Weighing sling MZ111600-002  |
| Batteries:   | 2 x CR2032 long life batteries   |