



The gently contoured baby tray gives generous as well as secure space for large and small infants. The fine graduation of 5 g all through the weighing range of 20 kg guarantees together with the auto-hold function precise and fast weighing results. The large white backlight LCD display makes the reading easy also in dimmed surroundings. Simply to operate by four clearly marked soft touch buttons. The fully integrated digital length measuring device will take precise readings of the infant up to 80cm by 2mm.

Technical specification:

Capacity	20 kg
Graduation	5 g
Length measuring range	400-800mm by 2mm
Functions	Tare, auto-hold, kg /lbs, cm/ Inch reading
Dimensions	600 x 280 x 135 mm
Weight approx.	3,5 kg
Power supply	4 x AA batteries or mains operated
Mains adapter	Optionally

