

## Column and floor scales

0.0 kg



The perfect scale for daily use at Hospitals, Clinics and GPs offices. Strong and sturdy construction with a total capacity of 250 kg and a fine graduation of 50 g. The turning display allows the reading from both sides of the display. Mains (optionally) and battery operated for maximum flexibility. A BMI function is integrated as standard function.

### Technical specification:

Capacity: 250 kg  
Graduation: 50g<50kg>100 g  
Dimensions: 400 x 350 x 44 mm  
Weight: 3,8 kg

Function: BMI calculation, mother and child weighing, turning display, hold & auto-hold, auto-switch off, tare function. Main all off switch. Power supply: 4 x AA batteries life time approx. 4000 cycles. Mains optionally, levelling feet, tap on start function