



Electronic wheelchair scale with integrated ramps on both sides. Handles and two castors allow an easy handling of the scale. The all steel construction gives extra stability and takes weights up to 300kg. The indicator with blue backlit LCD and be fix to the wall with the supplied wall mounting bracket. Powered by mains and rechargeable batteries.

#### Technical specification:

Capacity:	300 kg
Graduation:	100 g
Dimensions:	1070 x 1050 x 70 mm
Weight;	45 kg
Platform size:	800 x 800 mm
Functions:	Automatic zero setting, tare, pre-tare, low battery indication, auto witch off, hold, BMI, RS232-Interface, kg/lbs selectable