

# Automatic dial baby and toddler scale

## M108800

**ADE**



Quick to use baby and toddler weighing scale with a clear weight indication by two red revolving indicator needles. Robust and simple to handle with a full metal housing and weighing mechanism for general use as well as at rural situations for field workers. With a removable baby weighing tray to weigh toddlers and children standing upright up to 20kg.

### Technical specification:

Capacity:	20 kg
Graduation:	50 g < 10 kg > 100 g
Functions:	Manual zero setting
Dimensions weighing tray	550 x 265 x 80 mm
Overall dimensions:	568 x 259 x 343 mm
Weight approx.:	4,8 kg