

Column and floor scales

0.0 kg



Practical and easy to used electronic floor scale for routine weight checks of patients, athletes or for any person concerned about their weight. The good-sized weighing surface will take small or large people giving them a very secure standing due to the tilt-proof scale body. The clear red dot matrix display allows the reading under all lighting conditions with a minimum of power consumption. The step on start together with an auto-hold function will give quick and consistent readings with a fine graduation of just 50 g. A none-slip robust rubber surface will add to the safe and secure sense when standing on the scale. At a weight of just under 1.6 kg the scale can easily be carried.

Technical specification:

Capacity:	200 kg /440 lbs
Graduation:	50 g / 0,1 lbs
Dimensions :	338 x 36 x 340 mm
Weight :	1,6 kg

Start via 'Step-on' or 'Tap-on' function, auto-hold function, automatic zero setting, auto switch off, slip-proof surface, kg/lbs switch-over. Power supply 4 x AAA batteries life time approx. 4000 cycles.