

Electronic floor scale with mother and child function M321600



Column and floor scales

0.0 kg

The fully automatic mother and child weighing function will allow determining the weight of new born and small children fast and effortless in the arms of a nurse, parent or field worker just by starting the scale via the tap on function. A lightweight splash proof carry bag for transport is supplied as standard. A fully sealed on and off switch is placed beneath the scale to avoid unwanted battery use in transport. The light-coloured slip resistant platform protects the surface from unwanted heating up in bright sun light. The default setting will start the scale in the mother and child weighing mode. A manual switch over to standard weighing is quickly possible by pressing just one button.



Function: Mother and child weighing, automatic zero setting, auto-hold function, auto switch off, levelling feet, slip-proof surface

Power supply: 4 x AA batteries life time approx. 5000 cycles.

Tap on start function.

Battery operation: 180sec. / Mains power 7 Minutes

Auto-switch off & Main all off switch

Optional: ADE international mains power adapter H2870-006

Technical specification:

Capacity:	250 kg
Graduation:	100 g
Accuracy:	+/- 100g
Dimensions:	346 x 398 x 40-60 mm – levelling feet
Weight:	3,8 kg

