



This spacious platform scale has room for larger or obese patients and takes weights up to 300 kg with an accuracy of 100 g. The integrated BMI function allows quick results on the nutritional standard of the patient. The handle allows easy transport.

Technical specification:

Capacity:	300 kg
Graduation:	100 g
Functions:	Tare & pre-tare, hold, automatic zero setting, low battery indication, auto switch off, BMI calculation, kg/lbs switchable, RS-232 interface
Dimensions:	520 x 520 x 55 mm
Weight:	18 kg
Platform size:	520 x 520 mm