

# Electronic chair weighing scale

M400660 NEW

**ADE**



Comfortable and practical to use chair scale. Swivelling foot and arm rests allow a comfortable access in and out of the chair. 4 steering wheels make the scale easily manoeuvrable and safe due to two brakes at the back wheels. The tare and a hold function help to handle all patients. Powered by rechargeable batteries and mains for the use at any place in Hospitals.

#### Technical specification:

Capacity:	250 kg
Graduation:	100 g
Dimensions:	570 x 1055 x 935 mm
Weight;	28 kg
Seat size:	420 x 480 mm / seat height 530 mm
Functions:	Automatic zero setting, tare, pre-tare, low battery indication, auto witch off, hold, BMI, RS232-Interface, kg/lbs selectable

