


## Special purpose scale PW501 Compact and nappy scale

**ADE**



Diaper and nappy scale with hygienic stainless-steel tray. The scale measures accurately to 1 g and allows to keep track of babies output. The tare function allows subtracting the weight of small containers. With only 570 g in weight it allows easy transportation and handling.

- LCD screen with 10mm digits
- Automatic switch-off
- Removable stainless-steel platform
- Weighing & taring function
- Automatic zero-setting

### Technical specification:

Capacity:	5 kg
Graduation:	1 g
Dimensions:	177x223x38 mm
Weight:	570 g
Platform size:	177x140 mm
Power supply:	4 x AAA batteries or mains optionally