

Electronic bed weighing scale

M600620

ADE

Bariatric, bed and lifter scales

0.0 kg



For the use of bedridden patient in intensive care units and an essential must in renal kidney units during dialysis treatment. The bed can be easily pushed on the four weighing shoe with the help of the integrated ramps. The instrument trolley allows comfortable storage of the weighing shoe and a safe transport.

Technical specification:

Capacity:	500kg total capacity patient weight 300kg
Graduation:	100 g
Dimensions:	600 x 550 x 920 mm
Weight:	33 kg
Power supply:	By mains and rechargeable batteries
Function:	Automatic zero setting, tare, pre-tare, low battery indication, auto witch off, hold, BMI, RS232-Interface

