



Professional veterinarian scale for small and large pets with a fine graduation of 50 g. The stainless-steel platform guarantees best hygienic standards. A hold function will optimize the reading of the results on the display. Operated by mains or 4 x AAA batteries all included.

#### Technical specification:

Capacity:	150 kg
Graduation:	50 g
Dimensions:	900 x 550 x 75 mm
Weight:	20 kg
Functions:	Automatic zero setting, tare, low battery indication, auto witch off, hold function, kg/lbs reading
Power supply:	Mains adapter or 4 x AAA batteries