

Baby length measuring board

MZ10040

ADE



Perfectly suited for the mobile and stationary use. The two pieces of the measuring board are quickly clipped together and again disassembled allowing a space-saving storages or transport. The foot stop and the measuring boards are hold firmly together by just one patent clip mechanism. A skin-friendly even surface allows an easy cleaning of all parts. The smoothly running foot stop glider, guided by a central rail measures accurately to 1 mm. With a measuring range up to 100 cm even larger infants can be measures. Manufactured out of unbreakable ABS material for long-lasting trouble-free use under heavy duty conditions. Optionally available lightweight carry bag ADE MZ10043

Technical specification:

Measuring range	100-1000 mm
Graduation	1 mm
Dimensions	1260 x 280 x 90 mm
Weight approx.	585 g

