ADE Germany GmbH – Neuer Höltigbaum 15 – 22143 Hamburg - Germany - www.ade-germany.de









Durable padded light weight carrying bag made of water-repellent nylon material to protect, keep and transport our international and class III baby weighing scales. Right option for the baby scales: ADE M1001000-01, M118000, M118600, M112600, M112600U and M119600. Provided with a length adjustable removable shoulder strap and a handgrip for conformable transport. The strong and quality zipper will tightly close the bag against dust.

Technical specifications

Size: 600 x 380 x 150 mm Weight: 600 g Capacity: 5 kg Color: blue Material: Polyester and nylon fabric