


Column and floor scales



Low level platform scale with a stable handrail for additional stability and comfort for the patient. The scale can take weights up to 300 kg and offers support to patients with handicaps. Two castors allow an easy transport of the unit. Powered by mains and rechargeable batteries.

Technical specification:

Capacity:	300 kg
Graduation:	100 g
Functions:	Tare & pre-tare, hold, automatic zero setting, low battery indication, auto switch off, BMI calculation, kg/lbs switchable, RS-232 interface
Dimensions:	520 x 520 x 1245 mm
Weight:	18,5 kg
Platform size:	520 x 520 mm