



The ideal scale for your pet at home or at the veterinarian office for quick and effortless weight control. The multiplex layered wooden platform is covered with a non-slip PVC lining for easy cleaning, safe and comfortable use. With its weight of just over 4 kg the scale can be effortlessly carried. An easy-to-use fast reacting tap-on function allows starting the scale by foot or via a slight touch while holding a pet in the arms. The large blue backlite LCD is clear and fast to read even under dimmed light conditions at a minimum power consumption. A quickly responding auto-hold function will help to determine the weight, which can be read off after the pet has left the scale. After use, the scale will switch off automatically.

#### Technical specification:

Capacity:	100 kg
Graduation:	100 g
Dimensions:	800 x 600 x 33 mm
Weight:	4,3 kg
Functions:	Automatic zero setting, low battery indication, tap-on-start, auto-hold and auto switch off
Power supply:	2 x AAA alkaline batteries

