

PRODUCT INFORMATION



Electronic Table Scales



Compact Scale – Series HW945

- Perfect scale for control weighing
- Large LCD Display with backlight, height of digits 25mm
- Easy-to-clean membrane keypad
- Strong plastic housing
- Stainless steel weighing platform
- Protection cover for membrane keypad and housing
- Tare, Multitare
- Minus weighing mode to determine the weight of the amount taken
- Counting function
- Mains and rechargeable battery, 230V AC and rechargeable battery with 100 hours working hours

Technical Specifications		
Model	HW945-3	HW945-30
Art. No.	2306	2307
Capacity	3 kg	30 kg
Graduation	0.1 g	1 g
Dimensions (WxDxH)	260 x 310 x 125 mm	
Platform	260 x 215 mm	
Dead Weight	3.6 kg	
Power supply	230 V and rechargeable battery	

Subject to technical modifications due to further development.