

PRODUCT INFORMATION



Mechanical Round Head Scales



Floor scale – Series R (approved)

- Big round head
- Housing and round head lacquered
- Weighing platform hammer tone lacquered
- Manual zero setting
- Manual zero setting
- Easy to read weighing results
- Repeated rotation of the weighing locator for reliable measuring up to the maximum capacity

Technical Specifications		
Model	R60	R150
Art. No.	2292-60	2292-150
Capacity	60 kg	150 kg
Graduation	20 g	50 g
Weighing platform (WxD)	600 x 650 mm	
Dimensions (WxDxH)	600 x 890 x 1240 mm	
Dead weight	95 kg	

Technical modifications might be due to technical improvements.