



Precision scale – series PFA

- LCD-Display with backlight, height of digits 17mm
- Easy-to-clean membrane keypad
- Strong plastic housing
- Stainless steel weighing platform
- Weighing, taring, counting, percent
- Freely programmable plus-minus check weighing
- Auto-Off-Function
- With draft shield (only model PFA1000)
- Mains and battery (6x 1,5V AA)

Technical Specifications

Model	Art. no.	Capacity	Graduation	Weighing platform	Dimensions (WxDxH)	Dead weight
PFA-1000	1881	1,000 g	0.01 g	Ø 130 mm	190 x 260 x 78 mm	1.7 kg
PFA-3000	1882	3,000 g	0.1 g	160x160 mm		1.5 kg
PFA-6000	1883	6,000 g				

Subject to technical modifications due to further development.