

Electronic bed weighing scale
M601020 CE class III approved

ADE



For the use of bedridden patients in intensive care units and an essential must in renal kidney units during dialysis treatment. The bed can be easily pushed on the two independent weighbridges with integrated ramps for effortless placement of beds. The inbuilt rechargeable battery will store the weight even if an unwanted power failure occurs.

Technical specification:

Capacity: 500kg total capacity patient weight 300kg
Graduation: 100 g
Dimensions: 600 x 550 x 920 mm
Weight: 33 kg
Functions: Automatic zero setting, tare, pre-tare and memory,
low battery indication, auto switch off, hold, BMI, RS232-Interface,
including mobile stand