

Electronic baby scale

M118000 CE class III approved

ADE



The gently contoured baby tray offers generous as well as secure space for small infants and even for seated toddlers. The fine graduation of 5 g up to 10 kg, together with the auto-hold function and an electronic dampening system guarantees precise and fast weighing results. The clear backlit display ensures easy reading also in dimmed surroundings. Simple to operate by four clearly marked soft touch buttons. Optionally a carrying bag MZ10060 available.

Technical specification:

Capacity:	20 kg
Graduation:	5g<10kg>10g
Functions:	Tare, hold, automatic zero setting, low battery indication, automatic switch off
Dimensions:	600 x 350 x 135 mm
Dimensions tray:	600 x 280 x 90 mm
Weight approx.:	3,2 kg
Power supply:	4 x AA batteries / mains optionally

