

Electronic baby scale

M118000-01 CE class III approved

ADE



The gently contoured baby tray offers generous as well as secure space for small infants and even for seated toddlers. The fine graduation of 5 g up to 10 kg, together with the auto-hold function and an electronic dampening system guarantees precise and fast weighing results. The clear backlit display ensures easy reading also in dimmed surroundings. Simple to operate by four clearly marked soft touch buttons. The optionally integrated digital length measuring device will take precise readings of infants up to 80 cm in 2 mm steps. The infant can be freely positioned inside the weighing tray, no additional handling necessary. The measuring results will be shown on the LCD display of the baby weighing scale.

Technical specification:

Capacity:	20 kg
Graduation:	5g<10kg>10g
Length measure:	400 – 800 mm by 2 mm graduation
Functions:	Tare, hold, automatic zero setting, low battery indication, automatic switch off
Dimensions:	600 x 350 x 135 mm
Dimensions tray:	600 x 280 x 90 mm
Weight approx.:	3,2 kg
Power supply:	4 x AA batteries / mains optionally

