



Wide body floor scale to take into account a growing population of increasing size and weight. Perfect for hospital and clinic use with its integrated BMI calculation function. The display reading direction can be turned to the staff at the touch of a button or towards the patient. The tare and hold function completes the professional weighing scale. The mother and child function allows weighing a child in the arms of a parent.

Technische Spezifikationen:

Capacity:	250 kg	Function: BMI calculation, mother and child weighing, turning display,
Graduation:	50g<50kg>100 g	hold & auto-hold, auto-switch off, tap on start function, main all off switch,
Dimensions:	380 x 305 x 44 mm	levelling feet
Weight:	4 kg	Power supply: 4 x AA batteries life time approx. 4000 cycles and mains
Optionally	Carrying bag MZ10062	CE certified after directive 2014/31/EU and 93/42/EEC