

Electronic wheel chair and chair scales

M501020 Electronic wheel chair scale

ADE



Low level platform scale with two ramps for convenient access for wheel chairs. The all steel construction provides extra stability and takes weight up to 300 kg. A cable remote display shows the weight clearly and can be mounted to the wall. Two castors and two handgrips permit a quick transport. The pre tare function will allow you to deduct a memorized weight just by pressing one button. Such the net weight of a patient can be determined. A hold and tare function will help to get a quick and precise readings of the patient's weight. Powered by mains and rechargeable batteries all included.

Technical specification:

Capacity:	300 kg
Graduation:	100 g
Dimensions:	800 x 800 x 44 mm
Weight:	45 kg
Platform size:	800 x 800 mm
Functions:	Automatic zero setting, tare, pre-tare, low battery indication, auto witch off, hold, BMI calculation