

# Electronic chair weighing scale with 4 steering wheels M400020-01 CE class III approved

**ADE**



ADE M400020

Comfortable and practical to use the chair scale. Swivelling foot and arm rests allow a comfortable access in and out of the chair. The four steering wheels makes the scale easily manoeuvrable and safe due to brakes at the back and front wheels. Tare and the hold function help to handle all patients. Powered by rechargeable batteries or mains for the use at any place.

## Technical specification:

Capacity:	250 kg
Graduation:	100 g
Seat size:	420 x 480 mm / seat height 530 mm
Functions:	Automatic zero setting, tare, pre-tare and memory, low battery, indication, auto switch off, hold, BMI, RS232-Interface



ADE M400020-01