

# Electronic baby and toddler scale

M112800

**ADE**



Comfortable to use baby and toddler weighing scale for the check-up of babies' weight at the GPs office or at home. Quickly convertible to a floor scale with a capacity of 50kg for toddlers and young children. The last weight memory function will help to calculate the intake of milk by the infant and helps to keep track of the weight development of new born babies. Perfectly suitable for the rent-out service at pharmacies. A carrying bag is optionally available.

#### Technical specification:

Capacity:	up to 50 kg
Graduation:	10 g
Dimensions:	630 x 350 x 125 mm
Weight approx.:	1,6 kg
Optional extra:	Carrying bag MZ10060
Functions:	Tare, autohold, automatic zero setting, units switch kg/lbs, low battery indication
Special function:	BMIF-last weight memory function

