

## Column and floor scales

0.0 kg



Electronic handrail scale with two side rails for additional support of patients with disabilities and problems to stand allow on a weighing scale. The low level platform allows an easy access. With the pre-tare function any additional weight can be automatically deducted such as a small stool. Powered by mains and rechargeable batteries. Two castors are fixed to the platform for mobility.

### Technical specification:

Capacity:	300 kg
Graduation:	100 g
Functions:	Tare & pre-tare, hold, automatic zero setting, low battery indication, auto switch off, BMI calculation, kg/lbs switchable, RS-232 interface
Dimensions:	520 x 520 x 1135 mm
Weight:	18,5 kg
Platform size	520 x 520 mm