

Special purpose scale  
MP901600 electronic organ scale

**ADE**



With an accuracy of 1g this completely sealed scale is designed for specific needs in hospitals and medical institutions. The stainless steel bowl completes this unit. Powered by inbuilt rechargeable batteries for mains independent use.

**Technical specification:**

Capacity:	30 kg	Function:	Tare, automatic zero setting, low battery indication,
Graduation:	1 g		back light LCD, auto switch off
Dimensions:	535 x 330 mm	Power supply:	mains and rechargeable batteries
Weight:	4,1 kg	Levelling feet	