

Electronic wheel chair scale M500660

ADE



Low level platform scale with one ramp for convenient access. The cable remote LCD display with blue backlit shows the weight clearly and can be mounted to the wall with the supplied support bracket. Two castors and one handle permit a quick transport. Powered by mains and rechargeable batteries, allowing the use of the scale at any place. A second ramp can be ordered as optional accessory MZ30010.

Technical specification:

Capacity:	300 kg
Graduation:	100 g
Dimensions:	915 x 800 x 44 mm
Weight:	22 kg
Platform size:	800 x 800 mm
Functions:	Automatic zero setting, tare, pre-tare, low battery indication, auto witch off, hold, BMI, RS232-Interface, kg/lbs selectable
Option:	M500660-02 with two side boundaries



Optional M500660-02