



The perfect choice for small and medium sized pets. The robust and easy to clean non-slip weighing platform has a size of 71 x 41 cm and offers well enough space for cats and dogs as well as other pets up to a total capacity of 100 kg. Our MV303600 is a practical help when monitoring the weight of your pet. The tare and hold function will assist you to get fast and accurate weighing results. The tare function will allow you to remove unwanted weight from the reading, like an additional cover. The Hold function will help to get fast a stable weight reading even if the animal is moving on the platform. The scale is used by veterinary practices, breeders and private pet owners at home.

Technical specification:

Capacity:	100 kg
Graduation:	50 g
Dimensions:	716 mm x 40 mm x 412 mm
Weight:	8 kg
Functions:	Automatic zero setting, low battery indication, auto-switch off, HOLD function, tare function, unit switch
Power supply:	2 AAA batteries