


Baby length measuring mat

MZ10027-1



Ideally suitable for the portable use to measure the length of babies and toddlers fast and comfortable in a lying down position. A fix head piece and a sliding foot stop will give accurate and clear readings. Made of skin-friendly, easy to clean hygienic material. Resistant against common soft disinfection fluids. Foldable for space saving transport. Reading in cm only.



Technical specification:

Measuring range	100-990mm
Graduation	5mm
Dimensions	1200 x 300 x 120 mm
Weight	approx. 516 g