



The new ADE MZ10042 portable stadiometer is ideally suitable for field workers as well as for PGs and hospital use. The height rod can be quickly taken apart into 4 pieces which can be stored within the platform. The two wall distance pieces keep the height rod steady and vertical at the wall without any additional fixing. The scale can be read from both sides at the clear red indicator point.

Technical specification:

Measuring range	15-210 cm
Graduation	1 mm
Dimensions	34,8 x 214 x 41 cm
Weight approx.	2,3 kg
Optional accessory	ADE carrying bag MZ10043

